Medtronic CoreValve® System:

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
Elements Below Reflect Indications For Use According To The CE Mark

Diagnostic Findings	Non-Invasive		Angiography				Selection Criteria	
	Echo	CT/MRI	LV	Ao Root	CAG	Vascular	Recommended	Not Recommended
Atrial or Ventricular Thrombus	X						Not Present	Present
Sub Aortic Stenosis	X	X	X				Not Present	Present
LV Ejection Fraction	X		X				≥ 20%	< 20% Without Contractile Reserve
Mitral Regurgitation	X						≤ Grade 2	> Grade 2 Organic Reason
Vascular Access Diameter		X				X	≥ 6 mm Diameter	< 6 mm Diameter
Aortic and Vascular Disease		X				X	None to Moderate	Severe Vascular Disease
Indications For 26 mm CoreValve Device								
Annulus Diameter	X	X					20-23 mm	< 20 mm or > 23 mm
Ascending Aorta Diameter		X		X			≤ 40 mm	> 40 mm
Indications For 29 mm CoreValve Device								
Annulus Diameter	X	X					24-27 mm	< 24 mm or > 27 mm
Ascending Aorta Diameter		X		X			≤ 43 mm	> 43 mm

- [A] Annulus Diameter
- [B] Sinus of Valsalva Width
- [C] Ascending Aorta Diameter
- [D] Sinus of Valsalva Height
- [E] Frame Height (≈ 5 cm)

General Medical Guidance For Use Of CoreValve*

Diagnostic Findings	Non-Invasive		Angiography				Selection Criteria	
	Echo	CT/MRI	LV	Ao Root	CAG	Vascular	Recommended	Moderate-High Risk
LV Hypertrophy	X	X					Normal to Moderate 0.6 - 1.6 cm	Severe ≥ 1.7 cm
Coronary Artery Disease		X			X		None, Mid or Distal >70%	Proximal Lesions > 70%
Aortic Arch Angulation		X				X	Large Radial Turn	Sharp Turn
Aortic Root Angulation		X				X	< 30 Degrees	30 - 45 Degrees
Aortic and Vascular Disease		X				X	No or Light Vascular Disease	Moderate Vascular Disease
Vascular Access Diameter		X				X	> 6 mm	Calcified and Tortuous < 7 mm
Anatomic Considerations For 26 mm CoreValve Device								
Sinus of Valsalva Width	X	X		X			≥ 27 mm	< 27 mm
Sinus of Valsalva Height	X	X		X			≥ 15 mm	< 15 mm
Anatomic Considerations For 29 mm CoreValve Device								
Sinus of Valsalva Width	X	X		X			≥ 29 mm	< 29 mm
Sinus of Valsalva Height	X	X		X			≥ 15 mm	< 15 mm

* General medical guidance reflects the experience to date with the product, but final judgment remains with the implanting physician(s).

Consult with a certified proctor to determine if your patient is Moderate-High Risk.

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