

# AZURE™ S SR MRI SURESCAN™

Model W3SR01

## Physical characteristics

### Physical characteristics

Volume <sup>a</sup>	12.25 cm <sup>3</sup>
Mass	22.5 g
H x W x D <sup>b</sup>	42.6 mm x 50.8 mm x 7.4 mm
Radiopaque ID <sup>c</sup>	RNA
Surface area of titanium device can	33.48 cm <sup>2</sup>
Materials in contact with human tissue <sup>d</sup>	Titanium, polyurethane, silicone rubber
Battery	Lithium-hybrid CFx silver vanadium oxide

<sup>a</sup> Volume with connector holes unplugged.

<sup>b</sup> Grommets may protrude slightly beyond the can surface.

<sup>c</sup> The radiopaque ID, which includes a Medtronic-identifier symbol, can be viewed in a fluoroscopic image of the device.

<sup>d</sup> These materials have been successfully tested for the ability to avoid biological incompatibility. The device does not produce an injurious temperature in the surrounding tissue during normal operation.

### Replacement indicators

Recommended Replacement Time (RRT)	≤ 2.63 V on 3 consecutive daily automatic measurements
Elective Replacement Indicator (ERI)	3 months after RRT
End of Service (EOS)	3 months after ERI

## Pacing parameters

### Modes, rates, and intervals

Parameter	Programmable values
Mode	VVIR ⚡; VVI; VOO; OVO
Lower Rate <sup>a</sup>	30; 35 ... 60 ⚡; 70; 75 ... 150 min <sup>-1</sup> (±2 min <sup>-1</sup> )

<sup>a</sup> The corresponding Lower Rate Interval can be calculated as follows:  
Lower Rate Interval (ms) = 60,000/Lower Rate.

### RV parameters

Parameter	Programmable values
RV Amplitude	0.5; 0.75 ... 1.25 V (+0.125 V/-33%) 1.50 ... 3.5 ⚡ ... 5; 5.5; 6; 8 V (+15%/-33%) <sup>a</sup>
RV Pulse Width	0.03; 0.06 ms (± 10 μs); 0.1; 0.2; 0.3; 0.4 ⚡ ... 1.5 ms (± 25 μs)



- Bluetooth® wireless telemetry
- Approved for 1.5T and 3T MRI use

Medtronic

## RV parameters, cont'd.

Parameter	Programmable values
RV Sensitivity <sup>b</sup>	0.45; 0.60 mV ( $\pm 50\%$ ); 0.90 $\diamond$ ; 1.20; 2.00; 2.80; 4.00; 5.60; 8.00; 11.30 mV ( $\pm 30\%$ ) Unipolar: 2.80 $\diamond$ mV Bipolar: 0.90 $\diamond$ mV
RV Pace Polarity	Bipolar; Unipolar
RV Sense Polarity	Bipolar; Unipolar
RV Lead Monitor	Monitor Only; Adaptive
Min Limit	200 $\diamond$ ; 300; 400; 500 $\Omega$
Max Limit	1,000; 1,500; 2,000; 3,000 $\diamond \Omega$

<sup>a</sup> When RV Amplitude is 8 V, RV Pulse Width must be less than 1.3 ms.

<sup>b</sup> This setting applies to all sensing in this chamber for both tachyarrhythmia detection and bradycardia pacing operations.

## RV Capture Management™ parameters

Parameter	Programmable values
RV Capture Management™	Adaptive $\diamond$ ; Monitor; Off
RV Amplitude Safety Margin	1.5x; 2.0x $\diamond$ ; 2.5x; 3.0x
RV Minimum Adapted Amplitude	1.0; 1.5; 2.0 $\diamond$ ; 2.5; 3.0; 3.5 V
RV Acute Phase Remaining	Off; 30; 60; 90; 120 $\diamond$ ; 150 days

## Blanking periods

Parameter	Programmable values
V. Blank Post VP	150; 160 ... 200 $\diamond$ ... 320 ms (+5; -30 ms)
V. Blank Post VS	120 $\diamond$ ; 130 ... 170; 200; 220; 250; 280; 300; 320 ms (+2; -30 ms)

## Rate response pacing parameters

Parameter	Programmable values
Upper Sensor Rate	80; 85 ... 130 $\diamond$ ... 175 min <sup>-1</sup> ( $\pm 2$ min <sup>-1</sup> )
ADL Rate	60; 65 ... 95 $\diamond$ ... 170 min <sup>-1</sup> ( $\pm 2$ min <sup>-1</sup> )
Rate Profile Optimisation	On $\diamond$ ; Off
ADL Response	1; 2; 3 $\diamond$ ; 4; 5
Exertion Response	1; 2; 3 $\diamond$ ; 4; 5
Activity Threshold	Low; Medium Low $\diamond$ ; Medium High; High
Activity Acceleration	15; 30 $\diamond$ ; 60 s
Activity Deceleration	Exercise $\diamond$ ; 2.5; 5; 10 min
ADL Set Point	5; 6 ... 40; 42 ... 80
UR Set Point	15; 16 ... 40; 42 ... 80; 85 ... 180

## Sleep parameters

Parameter	Programmable values
Sleep	On; Off $\diamond$
Sleep Rate	30; 35 ... 50 $\diamond$ ; 55; 60; 70; 75 ... 100 min <sup>-1</sup>
Bed Time	00:00; 00:10 ... 22:00 $\diamond$ ... 23:50
Wake Time	00:00; 00:10 ... 07:00 $\diamond$ ... 23:50

## MRI SureScan™ parameters

Parameter	Programmable values
MRI SureScan™	On; Off
MRI Pacing Mode	VOO; OVO
MRI Pacing Rate	60; 70; 75; 80 ... 120 min <sup>-1</sup>

## Additional pacing features

Parameter	Programmable values
Rate Hysteresis	Off $\diamond$ ; 30; 40 ... 80 min <sup>-1</sup>

## Tachyarrhythmia parameters

### Tachyarrhythmia detection parameters

Parameter	Programmable values
VT Monitor	Monitor $\diamond$ ; Off
VT Monitor Interval (Rate) <sup>a</sup>	280; 290 ... 360 $\diamond$ ... 500 ms
RV Sensitivity <sup>b,c</sup>	0.45; 0.60 mV ( $\pm 50\%$ ); 0.90; 1.20; 2.00; 2.80; 4.00; 5.60; 8.00; 11.30 mV ( $\pm 30\%$ ) Bipolar: 0.90 $\diamond$ mV Unipolar: 2.80 $\diamond$ mV

<sup>a</sup> The measured intervals are truncated to a 10 ms multiple (for example, 457 ms becomes 450 ms). The device uses this truncated interval value when applying the programmed criteria and calculating interval averages.

<sup>b</sup> This setting applies to all sensing in this chamber for both tachyarrhythmia detection and bradycardia pacing operations.

<sup>c</sup> The device complies with the requirements of ISO 14708-2 when the sensitivity threshold is programmed to 2.0 mV or higher.

## Data collection parameters

### Data collection parameters

Parameter	Programmable values
EGM 1 Source	Can to RVring; RVtip to RVring $\diamond$ ; RVtip to Can
EGM 1 Range	$\pm 1$ ; $\pm 2$ ; $\pm 4$ ; $\pm 8$ $\diamond$ ; $\pm 12$ ; $\pm 16$ ; $\pm 32$ mV
EGM 2 Source	Can to RVring; RVtip to RVring; RVtip to Can $\diamond$
EGM 2 Range	$\pm 1$ ; $\pm 2$ ; $\pm 4$ ; $\pm 8$ $\diamond$ ; $\pm 12$ ; $\pm 16$ ; $\pm 32$ mV
EGM 3 Source	RVtip to RVring; Can to RVring $\diamond$ ; RVTip to Can
EGM 3 Range	$\pm 1$ ; $\pm 2$ ; $\pm 4$ ; $\pm 8$ $\diamond$ ; $\pm 12$ ; $\pm 16$ ; $\pm 32$ mV
Monitored	EGM1 and EGM2 $\diamond$ ; EGM1 and EGM3; EGM2 and EGM3
Pre-arrhythmia EGM	Off $\diamond$ ; On – 1 month; On – 3 months; On Continuous
Device Date/Time <sup>a</sup>	(select Time Zone)
Holter Telemetry	Off $\diamond$ ; 0.5; 1; 2; 4; 8; 16; 24; 36; 46 hr
Wireless Telemetry with Monitor	On; Off

<sup>a</sup> The times and dates stored in episode records and other data are determined by the Device Date/Time clock.

## Medtronic CareAlert™ parameters

### Clinical management alerts

Parameter	Programmable values
Monitored VT Episode Detected	Off ; On

### Lead/Device integrity alerts

Parameter	Programmable values
Low Battery Voltage RRT	On ; Off

### Lead Impedance Out of Range

Parameter	Programmable values
Lead Impedance	
RV Pacing Enable	On ; Off
RV Pacing Less than	200 ; 300; 400; 500 Ω
RV Pacing Greater than	1,000; 1,500; 2,000; 3,000  Ω
Capture Management™ High Threshold	
High Threshold	
RV Capture Enable <sup>a</sup>	Off ; On

<sup>a</sup> If programmed to On, alert notification is sent if RV capture management™ has measured high thresholds for 3 consecutive days.

## System test parameters

### System test parameters

Parameter	Selectable values
<b>Pacing Threshold Test parameters</b>	
Test Type	Amplitude; Pulse Width
Chamber	RV
Decrement after	2; 3 ... 15 pulses
RV Pace Polarity	Unipolar; Bipolar
Mode <sup>a</sup>	VVI; VOO
Lower Rate	30; 35 ... 60; 70; 75 ... 150 min <sup>-1</sup>
RV Amplitude	0.25; 0.5 ... 5; 5.5; 6; 8 V
RV Pulse Width	0.03; 0.06; 0.1; 0.2 ... 1.5 ms
V. Pace Blanking	150; 160 ... 320 ms
<b>Sensing Test parameters</b>	
Mode <sup>a</sup>	VVI; OVO
Lower Rate	30; 35 ... 60; 70; 75 ... 120 min <sup>-1</sup>

<sup>a</sup> The selectable values for this parameter depend on the programmed pacing mode.

## Longevity

### Projected service life in years

Pacing	Pre-arrhythmia EGM storage <sup>a</sup>	500 Ω pacing impedance		600 Ω pacing impedance		900 Ω pacing impedance	
		2.5 V	3.5 V	2.5 V	3.5 V	2.5 V	3.5 V
VVI, 0%	Off	18.3	18.3	18.3	18.3	18.3	18.3
	On	18.2	18.2	18.2	18.2	18.2	18.2
VVI, 15%	Off	17.5	16.7	17.6	16.9	17.8	17.3
	On	17.4	16.6	17.5	16.8	17.7	17.2
VVI, 50%	Off	15.8	13.9	16.1	14.4	16.7	15.4
	On	15.7	13.8	16.0	14.3	16.6	15.3
VVI, 100%	Off	13.9	11.2	14.4	11.9	15.4	13.3
	On	13.8	11.1	14.3	11.8	15.3	13.2

<sup>a</sup> The data provided for programming Pre-arrhythmia EGM storage to On are based on a 6-month period (two 3-month follow-up intervals) over the life of the device. Additional use of Pre-arrhythmia EGM storage reduces projected service life by approximately 13.6% or 1.6 months per year.

**Note:** These projections are based on typical shelf storage time (5 months). Assuming worst-case shelf storage time (18 months), longevity is reduced by approximately 7%.

The data are based on pacing outputs programmed to the specified amplitude and 0.4 ms pulse width and 60 min<sup>-1</sup> pacing rate. The service life of the device is affected by the programmed settings for certain features, such as Pre-arrhythmia EGM storage. Projected service life estimates are based on accelerated battery discharge data and device modeling as specified. These values should not be interpreted as precise numbers.

#### **IMPORTANT REMINDER**

This information is intended only for users in markets where Medtronic products and therapies are approved or available for use as indicated within the respective product manuals. Content on specific Medtronic products and therapies is not intended for users in markets that do not have authorization for use.

#### **Brief Statement**

See the MRI SureScan™ technical manual before performing an MRI scan and the device manual for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative or consult the Medtronic website at [www.medtronic.com](http://www.medtronic.com).



[www.medtronic.com/manuals](http://www.medtronic.com/manuals)

Consult instructions for use on this website. Manuals can be viewed using a current version of any major Internet browser. For best results, use Adobe Acrobat® Reader with the browser.

Medtronic and the Medtronic logo are trademarks of Medtronic.

™ Third party brands are trademarks of their respective owners.

All other brands are trademarks of a Medtronic company.

# Medtronic

#### **Europe**

Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 0 21 802 70 00  
Fax: +41 0 21 802 79 00

#### **United Kingdom/Ireland**

Medtronic Limited  
Building 9  
Croxley Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 0 1923 212213  
Fax: +44 0 1923 241004

UC201706570bEE ©2017 Medtronic.  
All Rights Reserved. Printed in Europe.

[medtronic.eu](http://medtronic.eu)