






NON-NAVIGATION THERAPEUTIC CATHETERS

EZ STEER® DS Bi-Directional Catheter








European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35U35R	7F	D-F 	8 mm	4	1-7-4 mm	115 cm
35U5JR	7F	F-J 	8 mm	4	1-7-4 mm	115 cm
35U33R*	7F	D-D 	8 mm	4	1-7-4 mm	115 cm
35U55R*	7F	F-F 	8 mm	4	1-7-4 mm	115 cm
35UJJR*	7F	J-J 	8 mm	4	1-7-4 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE® Generator

CELSIUS® THERMOCOOL® Uni-Directional Catheter




European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35Q13R	7F	B 	3.5 mm	4	2-5-2 mm	115 cm
35Q33R	7F	D 	3.5 mm	4	2-5-2 mm	115 cm
35Q53R	7F	F 	3.5 mm	4	2-5-2 mm	115 cm
35Q73R*	7F	J 	3.5 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE® Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

* For availability, please contact your local Biosense Webster representative.

NON-NAVIGATION THERAPEUTIC CATHETERS

CELSIUS® RMT THERMOCOOL® Catheter



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit*	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
CR7TCSIRT**	8F	≥8.5F	3.5 mm	4	2-5-2 mm	130 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE® Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

CELSIUS® 4 mm Catheter Thermocouple



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35G07R	6F	A ●	4 mm	4	2-5-2 mm	92 cm
35G17R	6F	B ●	4 mm	4	2-5-2 mm	92 cm
35G37R	6F	D ●	4 mm	4	2-5-2 mm	92 cm
36N07R	7F	A ●	4 mm	4	2-5-2 mm	115 cm
36N17R	7F	B ●	4 mm	4	2-5-2 mm	115 cm
36N27R	7F	C ●	4 mm	4	2-5-2 mm	115 cm
36N37R	7F	D ●	4 mm	4	2-5-2 mm	115 cm
36N47R	7F	E ○	4 mm	4	2-5-2 mm	115 cm
36N57R	7F	F ●	4 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE® Generator







* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

** For availability, please contact your local Biosense Webster representative.

NON-NAVIGATION THERAPEUTIC CATHETERS

CELSIUS® 4 mm Catheter Thermistor







European Ordering No.	Outer Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36O07R	7F	A 	4 mm	4	2-5-2 mm	115 cm
36O17R	7F	B 	4 mm	4	2-5-2 mm	115 cm
36O27R	7F	C 	4 mm	4	2-5-2 mm	115 cm
36O37R	7F	D 	4 mm	4	2-5-2 mm	115 cm
36O47R	7F	E 	4 mm	4	2-5-2 mm	115 cm
36O57R	7F	F 	4 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE® Generator

CELSIUS® 4 mm Braided Tip Catheter




European Ordering No.	Outer Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
34D37R	7F	D 	4 mm	4	2-5-2 mm	115 cm
34D47R	7F	E 	4 mm	4	2-5-2 mm	115 cm
34D57R	7F	F 	4 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE® Generator

NON-NAVIGATION THERAPEUTIC CATHETERS

CELSIUS FLTR® 8 mm Uni-Directional Catheter



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D135501	7.5F	≥8F	F	8 mm	4	0.5-5-2 mm	115 cm
D135502	7.5F	≥8F	XR	8 mm	4	0.5-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE® Generator

CELSIUS FLTR® 8 mm Bi-Directional Catheter



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D135901	7.5F	≥8F	D-XR	8 mm	4	0.5-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE® Generator

* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

NON-NAVIGATION THERAPEUTIC CATHETERS

CELSIUS® DS Catheter



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35S17R	7F	≥8F	B ● ●	8 mm	4	1-6-2 mm	115 cm
35S37R	7F	≥8F	D ● ●	8 mm	4	1-6-2 mm	115 cm
35S57R	7F	≥8F	F ● ●	8 mm	4	1-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE® Generator

CELSIUS® RMT Catheter



European Ordering No.	Outer Diameter	Sheath Fit*	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
CR7TCS4RT	7F	≥8F	4 mm	4	2-5-2 mm	125 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE® Generator

* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.