

## DECANAV<sup>®</sup> Catheter



Pairs at 10 mm (2-8-2 mm)



Catheter handle

European Ordering No.	Diameter	Curve	Tip Size	Number of Electrodes	Spacing	Insertion Length
R7D282CT	7F	D ●	2 mm	10	2-8-2 mm	115 cm
R7F282CT	7F	F ●	2 mm	10	2-8-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3434CT	Catheter to CARTO <sup>®</sup> 3 System PIU