NON AUTO-ID DIAGNOSTIC CATHETERS

HALO® XP Tricuspid Mapping Catheter





European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
36J13R	7F	Tricuspid	20	(20)-2-8-2 mm	110 cm
36J16R	7F	Tricuspid	20	2-13-2 mm	110 cm
36J14R	7F	Tricuspid	20	2-18-2-8-2 mm	110 cm

Associated Cable:

Catalog No.:	Connects	
2 x 39F41R ¹	ECG Recording System CARTO® 3 System	
D128624	Catheter to CARTO® 3 System	

