

## LASSO® 2515 Variable Mapping Catheter



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35O2515R	7F	D ●	15-25 mm	10	8 mm	115 cm
35T2515R	7F	D ●	15-25 mm	20	2-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F41R <sup>1</sup> or 39F54R <sup>1</sup>	ECG Recording System CARTO® 3 System
D128624	Catheter to CARTO® 3 System

## LASSO® Mapping Catheter

European Ordering No.	Diameter	Curve	Loop Diameter	Number of Electrodes	Spacing	Insertion Length
35O16R	7F	D ●	12 mm	10	4 mm	115 cm
35O26R	7F	D ●	15 mm	10	4.5 mm	115 cm
35O36R	7F	D ●	20 mm	10	6 mm	115 cm
35O46R	7F	D ●	25 mm	10	8 mm	115 cm
35O66R*	7F	D ●	35 mm	15	7 mm	115 cm
35T26R	7F	D ●	15 mm	20	4.5 mm	115 cm
35T36R	7F	D ●	20 mm	20	6 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F41R <sup>1</sup> or 39F54R <sup>1</sup>	ECG Recording System CARTO® 3 System
D128624	Catheter to CARTO® 3 System

<sup>1</sup> The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.  
\* For availability, please contact your local Biosense Webster representative.