

NAVIGATION THERAPEUTIC CATHETERS

NAVISTAR® THERMOCOOL® Uni-Directional Catheter



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
34H17M	8F	≥8F	B ●	3.5 mm	4	2-5-2 mm	115 cm
34H27M	8F	≥8F	C ●	3.5 mm	4	2-5-2 mm	115 cm
34H37M	8F	≥8F	D ●	3.5 mm	4	2-5-2 mm	115 cm
34H57M	8F	≥8F	F ●	3.5 mm	4	2-5-2 mm	115 cm
34HJ7M	8F	≥8F	J ●	3.5 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO® 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE® Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

NAVISTAR® THERMOCOOL® Bi-Directional Catheter



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36H33M	8F	≥8F	D-D ●●	3.5 mm	4	2-5-2 mm	115 cm
36H35M	8F	≥8F	D-F ●●	3.5 mm	4	2-5-2 mm	115 cm
36H55M	8F	≥8F	F-F ●●	3.5 mm	4	2-5-2 mm	115 cm
36H5JM	8F	≥8F	F-J ●●	3.5 mm	4	2-5-2 mm	115 cm
36HJJM	8F	≥8F	J-J ●●	3.5 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO® 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE® Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.