

PARAHISIAN® Catheter



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36L00Q	7F	D ●	1 mm	8	2-10-2-2-2-2 mm	90 cm

Associated Cable:

Catalog No.:	Connects
39F44R	ECG Recording System / CARTO® XP

ECG Recording System

Associated Cable:

Catalog No.:	Connects
39F44R CY1210CT	ECG Recording System / CARTO® XP System CARTO® 3 System

1 The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.