

# NON-NAVIGATION THERAPEUTIC CATHETERS

## THERMOCOOL® SF Uni-Directional Catheter\*



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D131601	8F	≥8F	B ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131602	8F	≥8F	D ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131603	8F	≥8F	F ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131604	8F	≥8F	J ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing

## THERMOCOOL® SF Bi-Directional Catheter\*



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
BDI35BBRT	8F	≥8F	B-B ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35BDRT	8F	≥8F	B-D ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35BFRT	8F	≥8F	B-F ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35DDRT	8F	≥8F	D-D ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35DFRT	8F	≥8F	D-F ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35DJRT	8F	≥8F	D-J ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35FFRT	8F	≥8F	F-F ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35FJRT	8F	≥8F	F-J ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35JJRT	8F	≥8F	J-J ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing

\* THERMOCOOL® SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

\*\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.