

NAVIGATION THERAPEUTIC CATHETERS

THERMOCOOL® SF NAV Uni-Directional Catheter*



2-5-2 mm

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D131801	8F	≥8F	B ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131802	8F	≥8F	D ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131803	8F	≥8F	F ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131804	8F	≥8F	J ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO® 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE® Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

* THERMOCOOL® SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

** Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

NAVIGATION THERAPEUTIC CATHETERS

THERMOCOOL® SF NAV Bi-Directional Catheter*




European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
BNI35BBH	8F	≥8F	B-B ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35BDH	8F	≥8F	B-D ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35BFH	8F	≥8F	B-F ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35DDH	8F	≥8F	D-D ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35DFH	8F	≥8F	D-F ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35DJH	8F	≥8F	D-J ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35FFH	8F	≥8F	F-F ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35FJH	8F	≥8F	F-J ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35JJH	8F	≥8F	J-J ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO® 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE® Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

* THERMOCOOL® SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

** Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

NAVIGATION THERAPEUTIC CATHETERS

THERMOCOOL® SF NAV Uni-Directional Catheter* with curve visualization



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D131501	8F	≥8F	B ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131502	8F	≥8F	D ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131503	8F	≥8F	F ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131504	8F	≥8F	J ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3434CT	Catheter to CARTO® 3 System PIU
39E43R	PIU to Stockert EP Shuttle
D129704	PIU to SMARTABLATE® Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

THERMOCOOL® SF NAV Bi-Directional Catheter* with curve visualization



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
BNI35BDCT	8F	≥8F	B-D ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35BFCT	8F	≥8F	B-F ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35FJCT	8F	≥8F	F-J ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3434CT	Catheter to CARTO® 3 System PIU

* THERMOCOOL® SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

** Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.