

NAVIGATION THERAPEUTIC CATHETERS

THERMOCOOL SMARTTOUCH® SF Uni-Directional Catheter*

UNI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Spacing	Insertion Length
D134701	8F	≥8.5F	D ●	3.5 mm	1-6-2 mm	115 cm
D134702	8F	≥8.5F	F ●	3.5 mm	1-6-2 mm	115 cm
D134703	8F	≥8.5F	J ●	3.5 mm	1-6-2 mm	115 cm

THERMOCOOL SMARTTOUCH® SF Bi-Directional Catheter*

BI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Spacing	Insertion Length
D134801	8F	≥8.5F	D-D ●●	3.5 mm	1-6-2 mm	115 cm
D134802	8F	≥8.5F	F-F ●●	3.5 mm	1-6-2 mm	115 cm
D134803	8F	≥8.5F	J-J ●●	3.5 mm	1-6-2 mm	115 cm
D134804	8F	≥8.5F	F-J ●●	3.5 mm	1-6-2 mm	115 cm
D134805	8F	≥8.5F	D-F ●●	3.5 mm	1-6-2 mm	115 cm

THERMOCOOL SMARTTOUCH® Uni-Directional Catheter

UNI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D133601	8F	≥8.5F	D ●	3.5 mm	6	1-6-2 mm	115 cm
D133602	8F	≥8.5F	F ●	3.5 mm	6	1-6-2 mm	115 cm
D133603	8F	≥8.5F	J ●	3.5 mm	6	1-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3434CT	Catheter to CARTO® 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

* THERMOCOOL® SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

** Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

NAVIGATION THERAPEUTIC CATHETERS

THERMOCOOL SMARTTOUCH® Bi-Directional Catheter



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D132701	8F	≥8.5F	D-D ● ●	3.5 mm	6	1-6-2 mm	115 cm
D132702	8F	≥8.5F	F-F ● ●	3.5 mm	6	1-6-2 mm	115 cm
D132703	8F	≥8.5F	J-J ● ●	3.5 mm	6	1-6-2 mm	115 cm
D132704	8F	≥8.5F	F-J ● ●	3.5 mm	6	1-6-2 mm	115 cm
D132705	8F	≥8.5F	D-F ● ●	3.5 mm	6	1-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3434CT	Catheter to CARTO® 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator

Associated Tubing

Catalog No.:	Description
CFT003	COOLFLOW® Pump Tubing
SAT001	SMARTABLATE® System Irrigation Tubing Set

* THERMOCOOL® SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

** Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.