

## WEBSTER® CS Uni-Directional Catheter



2-8-2 mm



catheter handle

European Ordering No.	Diameter	Curve	Tip Electrodes	Number of Electrodes	Spacing	Insertion Length	Auto ID	Compression Coil
D135303	6F	CS-D	2 mm	10	2-8-2 mm	115 cm	Yes	Yes
D135304	6F	CS-F	2 mm	10	2-8-2 mm	115 cm	Yes	Yes

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO® 3 System PIU
CY1212CT	Catheter to CARTO® 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox



## WEBSTER® CS Bi-Directional Catheter



2-8-2 mm



catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
36G35QCT	7F	D-F 	10	2-8-2 mm	115 cm
36G5JMCT	7F	F-J 	10	2-8-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO® 3 System PIU
CY1212CT	Catheter to CARTO® 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

## WEBSTER® Quadrapolar Catheter Deflectable



European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
36B38RCT	6F	D ●	4	2-5-2 mm	92 cm
36B33RCT	6F	D ●	4	5 mm	92 cm
36C38RCT	7F	D ●	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO® 3 System PIU
CY1212CT	Catheter to CARTO® 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

## WEBSTER® Decapolar Catheter Deflectable



European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
35X33QCT	6F	D ●	10	2-8-2 mm	92 cm
36G33QCT	7F	D ●	10	2-8-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO® 3 System PIU
CY1212CT	Catheter to CARTO® 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

## WEBSTER® HIS Catheter



Catheter handle

European Ordering No.	Diameter	Curve	Poles	Spacing	Insertion Length
D1085413	5F	HIS-F	4	5 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO® 3 System PIU
CY1212CT	Catheter to CARTO® 3 System PIU
D128709	Auto ID Catheters up to 10 poles to the EP recording system pinbox

## WEBSTER® Quadrapolar Catheter Fixed



5 mm



Catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
37C53RCT	5F	F ●	4	5 mm	110 cm
37C03RCT	5F	A ●	4	5 mm	110 cm
F6QF005CT	6F	F ●	4	5 mm	115 cm
F6QA005CT	6F	A ●	4	5 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO® 3 System PIU
CY1212CT	Catheter to CARTO® 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

## WEBSTER® Decapolar Catheter Fixed



2-8-2 mm



Catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
37S94RCT	5F	P	10	2-8-2 mm	60 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO® 3 System PIU