

Coronary Balloons



Features

- Outstanding performance for kissing balloon under 6F guiding catheter.
- Patented hydrophilic coating on distal shaft and Hypotube proximal shaft provide outstanding deliverability.
- Three-directional memory folding of the balloon allows of size minimization with great withdrawal capability, and ensures the good performance of repeated dilatation.
- Soft tapered tip provides smooth transition from guidewire to balloon catheter to avoid "fish mouth", and minimize the vessel damage.
- Short shoulder design at both ends of the balloon ensures the precise dilatation, and reduces the damage to normal vessel next to the lesion.

Tadpole™

Coronary Dilatation Catheter



Technical Parameters

Catheter Design	Rapid Exchange
Balloon Compliance	Semi-compliant
Marker	2
Distal Shaft O.D.	2.7F
Proximal Shaft O.D.	1.9F
Tip Profile	0.017"
Minimum Guiding Catheter Compatibility	5F (0.056")
Guide Wire Compatibility	0.014"
Effective Length	142cm
Nominal Pressure	8 atm (φ2.25mm-3.0mm) / 6 atm (φ3.5mm-4.0mm)
Rated Burst Pressure	16 atm (φ2.25mm-3.75mm) / 14 atm (φ4.0mm)

Ordering Information

Balloon Diameter (mm)	Balloon Length (mm)									
	8	9	10	12	15	16	20	24	25	30
2.5	LPCR 25008	LPCR 25009	LPCR 25010	LPCR 25012	LPCR 25015	LPCR 25016	LPCR 25020	LPCR 25024	LPCR 25025	LPCR 25030
2.75	LPCR 27008	LPCR 27009	LPCR 27010	LPCR 27012	LPCR 27015	LPCR 27016	LPCR 27020	LPCR 27024	LPCR 27025	LPCR 27030
3.0	LPCR 30008	LPCR 30009	LPCR 30010	LPCR 30012	LPCR 30015	LPCR 30016	LPCR 30020	LPCR 30024	LPCR 30025	LPCR 30030
3.5	LPCR 35008	LPCR 35009	LPCR 35010	LPCR 35012	LPCR 35015	LPCR 35016	LPCR 35020	LPCR 35024	LPCR 35025	LPCR 35030
4.0	LPCR 40008	LPCR 40009	LPCR 40010	LPCR 40012	LPCR 40015	LPCR 40016	LPCR 40020	LPCR 40024	LPCR 40025	LPCR 40030

Notice: The sizes with * mark require a minimum quantity of 10 pcs per order. Please contact local sales for details before placing order.